

# R-1/R-2 FLOW CONTROL UNIT



The AEI Technologies R-1 and R-2 Flow Controllers are commonly used as an integral part of any gas measurement system. Either Flow Controller is frequently used in conjunction with the S-3A Oxygen analyzer and CD-3A Carbon Dioxide analyzer. The Model R-1 includes a flow meter, pump, and a needle valve in a single cabinet. Model R-2 (shown above) includes 2 flow meters, pump, and 2 needle valves.



**AEI Technologies, Inc.**  
520 East Ogden Avenue • Naperville, IL 60563  
[1] 800-793-7751 • [1] 630-548-3545 • f:[1] 630-548-3546  
[sales@aeitechnologies.com](mailto:sales@aeitechnologies.com) • [www.aeitechnologies.com](http://www.aeitechnologies.com)